



Spiration[®] Valve System (SVS) Patient Management Recommendations



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION




Patient Management Recommendations

Overview

Proper patient management is of the utmost importance for endobronchial valve treatment. Due to the complexity of this procedure, each step of the patient management process, from the initial patient screening to the post-procedure re-evaluation, requires a thorough assessment to ensure successful outcomes.

Patient Process

The information provided in this guide is based upon findings established by the EMPROVE clinical trial in combination with clinical subject matter expert's recommendations. This information is not meant to replace patient-specific clinical judgment. The following patient management workflow is recommended for optimal procedure results:



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Table of Contents

| | |
|--|----|
| ■ Initial Patient Screening | 4 |
| ■ Medical Patient Selection | 5 |
| ■ Radiographic Patient Selection | 6 |
| ■ Pre-Admission Protocol | 7 |
| ■ Procedure Day Protocol | 8 |
| ■ Post-Procedure Management | 11 |
| ■ Patient Follow-Up | 14 |
| ■ Post-Procedural Re-Evaluation | 16 |
| ■ References | 17 |

INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP


POST-PROCEDURAL
RE-EVALUATION



Initial Patient Screening

Medical History and Physical Exam

- Patient is 18 years of age or older
- Patient is not an active smoker
- Patient has GOLD stage III or IV COPD (emphysema)
- Patient has a BMI > 15 and < 35 kg/m²
- Patient is on optimized medical management
- Patient has had no prior, ipsilateral lung volume reduction surgery or major lung procedures (e.g. lobectomy, segmentectomy, transplant)
- **Pulmonary function testing**
 - Spirometry (Post-BD)
 - FEV₁ > 15 and $\leq 45\%$
 - If change in FEV₁ > 200 mL post-bronchodilator, patient not considered a good candidate for valve therapy
 - Body Plethysmography
 - RV $\geq 150\%$
 - TLC $\geq 100\%$
 - DLCO
 - Should be $> 20\%$



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP


POST-PROCEDURAL
RE-EVALUATION



Medical Patient Selection

Medical Work-Up

- Patient must perform a 6-minute walk distance of $\geq 100\text{m}$
- Patient must have severe dyspnea which is defined as a mMRC ≥ 2
- Patient meets the criteria of the ATS/ERS Guidelines for Management of Stable COPD
- Complete medical history to rule out patient has no co-existing major medical disease, alcoholism, or drug abuse. This includes neurological or musculoskeletal conditions that may interfere with patient selection testing
- Patient does not have clinically significant bronchitis
- Patient does not have an active asthma component to their disease or requires more than 20 mg of prednisone daily
- Patient has had a complete cardiac workup including evaluation for severe pulmonary hypertension based (Pa pressure $> 50\text{ mmHg}$)
- Patient is classified as ASA Class greater than P4 including presence of co-morbidity that could significantly increase the risk of a bronchoscopy procedure
- Assess patient for Alpha-1 antitrypsin deficiency
- **Arterial Blood Gas level**
 - Collected on room air, discontinue supplemental O_2 for 10 minutes prior to sampling if needed
 - $\text{PaCO}_2 < 55\text{ mm HG}$
 - $\text{PaO}_2 > 45\text{ mm HG}$



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Radiographic Patient Selection

Radiographic Assessment

- Presence of any lung nodule or other clinically significant radiographic findings have been ruled out or treated
- Verify that there is no giant bulla (> 30% volume in either lung) present
- Based on the EMPROVE clinical trial result the patient should meet the following criteria:
 - Patient should have severe emphysema and high heterogeneity (measured at -920 HU) defined as:
 - A target lobe with $\geq 40\%$ emphysema involvement
 - ≥ 10 percentage points disease severity difference with the ipsilateral lobe
 - Patients target lobe and ipsilateral lobe will be separated with an intact fissure
 - An intact fissure will be estimated visually to be ≥ 0 90% complete with no segmental vessels crossing from one lobe to the adjacent lobe after viewing the HRCT in 3 dimensions (use of the SeleCT report is the preferred method for determining fissure completeness)
 - Patient does not have a diffuse pattern of emphysema
- Perfusion Analysis (at physician's discretion)
 - Help to determine target lobe
 - Target the lobe with the lowest perfusion



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

**RADIOGRAPHIC
PATIENT SELECTION**

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP


POST-PROCEDURAL
RE-EVALUATION



Pre-Admission Protocol

Additional Tests Done within 30 Days of Procedure

- Verify H&P is current
- Verify smoking cessation > 4 months
- Required lab work
 - CBC
 - Chem 7
 - PT/PTT/INR (at physician's discretion)
 - If patient is on anticoagulant, it should be confirmed with the primary care physician or cardiologist that anticoagulants have been safely stopped 5-7 days before procedure
- Patient has had no hospitalizations for COPD exacerbation or respiratory infections in the past 3 months
- **Cardiopulmonary assessment**
 - Oxygenation and ventilation
 - Oxygen saturation at rest and exercise
 - Arterial blood gas
 - Use of CPAP or BiPAP
 - Cardiac function
 - Echocardiogram
 - No congestive heart failure (LVEF < 45%)
 - Pulmonary hypertension
 - sPAP > 50mmHg on right heart catheterization



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP


POST-PROCEDURAL
RE-EVALUATION



Procedure Day Protocol

Endobronchial Valve Treatment Preparations

- Prophylactic antibiotics and steroids (at physician's discretion)
 - Decreases incidence of AECOPD or Infection
 - Day of procedure and continue for 5 days
- Provide nebulized bronchodilators 20 minutes prior to case
 - Albuterol/Ipratropium bromide
- Consider a decrease in tidal volume by 20% after last valve placement to avoid non-targeted lobe overdistention until spontaneous ventilation resumes
- Provide post-procedural bronchodilators immediately after procedure in recovery area and then every 4–6 hours as needed
- **Anesthesia concerns**
 - ASA statement on NORA locations
 - Adequate space for all equipment and personnel
 - Reliable oxygen source and backup supply
 - Adequate and reliable suction source should meet OR standards
 - If inhaled anesthetics used, reliable scavenging system for waste anesthesia gases
 - Adequate illumination
 - Adequate electrical outlets with emergency backup power
 - Standard and emergency anesthesia supplies and drugs
 - Code cart and defibrillator
 - Back-up manual ventilation equipment



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Procedure Day Protocol (continued)

Endobronchial Valve Treatment Preparations

- Monitoring equipment that meets ASA standards
 - ECG continuously displayed
 - SpO₂ with audible tone with pitch determined by SpO₂
 - End-tidal carbon dioxide
 - BP at least every 5 minutes
 - Mechanical ventilator or anesthesia machine must have audible disconnect alarm
 - Temperature (ability to monitor if needed)
- General anesthesia (GA) vs. sedation, based on:
 - Both are acceptable
 - GA is more common and recommended
 - Determination of GA vs. sedation
 - Local preference and experience
 - Pre-anesthesia patient evaluation
- General Anesthesia
 - Total Intravenous Anesthesia (TIVA)
 - Does not require anesthesia machine
 - May be easier to implement in an endoscopy suite
 - Use of short acting medication can result in quick emergence with minimal hangover effect
 - If no anesthesia machine, a standard ventilator can be used
 - Requires one of more infusion pumps
 - Possible Anesthesia Options
 - Propofol
 - Opioid
 - Remi-fentanyl
 - Fentanyl
 - Sufentanil
 - Topical Anesthesia – 1% lidocaine
 - Others: Dexmedetomidine and ketamine

INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Procedure Day Protocol (continued)

Endobronchial Valve Treatment Preparations

- Inhalational anesthetic
 - Volatile agents cause bronchodilation
 - Reliable
 - Easy to titrate
 - Requires use of an anesthesia machine and scavenging of waste anesthesia gases
 - Insertion/removal of bronchoscope, airway suction
 - possible fluctuating levels of anesthesia
 - potential pollution of room
- Airway Management
 - Endotracheal tube
 - Reliable oxygenation and ventilation
 - Minimum size 8 mm preferable 8.5–9 mm (air leaks require a minimum of 8.5 mm)
- Neuromuscular Blockade
 - Used to facilitate tracheal intubation during induction of anesthesia
- Topical anesthesia
 - Will decrease cough and anesthetic requirements
- Anticholinergics may decrease secretions
 - Glycopyrrolate or atropine
 - May increase risk of tachyarrhythmias or urinary retention
- Sedation
 - Level of monitoring the same as GA
 - Topical anesthesia
 - Upper airway
 - Tracheobronchial tree
 - Technique
 - Benzodiazepines
 - Opioids (fentanyl, sufentanil, remifentanil)
 - Propofol
 - Dexmedetomidine
 - Ketamine
 - Plan for managing apnea/hypoventilation

INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Post-Procedure Management

Pneumothorax Management

- Monitoring for pneumothorax
 - Chest x-ray
 - Within 4 hours of procedure
 - Daily until discharge
 - Chest ultrasound
 - Optional, if local expertise present
 - Immediate post-procedure assessment for lobar occlusion and check for lung sliding. In absence of sliding, re-inspection before waking up the patient
 - Chest CT scan
 - Low threshold for CT imaging if pneumothorax is suspected and patient has rapid major lobar collapse to rule out associated PTX
- Pneumothorax cart/tray should always be in patient room after the procedure
 - Always be prepared to place a chest tube within 15-30 minutes during hospital stay
 - Needle/catheter for needle decompression (16Ga)
 - Seldinger chest tubes 14Fr
 - Surgical chest tubes and tray 20–28Fr for “open” chest tubes
 - Gloves, scalpel, suture, dressing materials, pleuro vac



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Post-Procedure Management (continued)

Post-Operative Day 0

- Hospital admission
 - Hospital stay is based upon physician’s discretion and should be considered if there is a persistent chest discomfort or pain on the treated side
- Place patient on telemetry monitoring and pulse oximetry
- Supplemental oxygen to maintain SpO₂ >92%
- Maintain bedrest
- Nebulized Ipratropium/Albuterol q4h, prn
- Restart home medications, LAMA/LABA/ICS
- DVT prophylaxis

Post-Operative Day 1–3

- Daily morning portable chest x-ray
- Day #1 ambulate in room with bathroom privileges
- Day #2-3 ambulate in hallway with PT assistance with telemetry
- Monitor for complications
 - Pneumothorax
 - COPD exacerbation symptoms
 - Valve settling
 - Rapid targeted lobe atelectasis
 - Consider obtaining a chest CT to rule out early minor pneumothorax
 - Increasing hypoxemia
 - Rule out pneumothorax
 - AECOPD
 - PNA
 - Airway kinking
 - CHF
 - Alveolar hemorrhage
 - Re-expansion edema of ipsilateral non-targeted lobe
 - Consider HFNC or NIV
 - Consider ICU care
- Increasing CO₂ retention
 - Rule out AECOPD
 - Increase BDA, systemic steroids
 - Consider BIPAP
 - ICU care

INITIAL PATIENT SCREENING

MEDICAL PATIENT SELECTION

RADIOGRAPHIC PATIENT SELECTION

PRE-ADMISSION PROTOCOL

PROCEDURE DAY PROTOCOL

POST-PROCEDURE MANAGEMENT

PATIENT FOLLOW-UP

POST-PROCEDURAL RE-EVALUATION



Post-Procedure Management (continued)

Day of Discharge

- PA and lateral chest x-ray
- Review discharge medications
- Give patient, a family member and/or spouse, a wallet card and emergency information
 - Number of valves placed and their location
 - Shows what lobe was treated
 - Include treating physicians name and contact number
 - Includes MRI conditional information
- Schedule patient to return for follow-up visit with treating pulmonary physician in a week
- Inform patient to call if he/she has any increased symptoms of chest pain, shortness of breath, cough, purulent sputum, fever, chills or hemoptysis
- Inform patient to go to the nearest ER if symptoms are severe or persist
- If patient is from a remote area or long distance from the hospital, consider having them stay with a family member or friend for the first few days after discharge to be close to treating facility.
 - If patient does not have friends or family in the area, consider having them stay in a local hotel

INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Patient Follow-Up

Short-Term

- Call patient daily for first 10 days post procedure
- Clinic follow-up in 1 week
 - Perform chest x-ray
 - 6MWD
- Clinic follow-up in 1 month
 - H&P
 - HRCT between 30–45 days to check the positioning of the valves and lobar collapse (at physician's discretion)
- Clinic follow-up in 3 months
 - Chest x-ray
 - PA/LAT
 - PFTS
 - 6MWD
 - SGRQ, mMRC, CAT and BODE



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Patient Follow-Up (continued)

Long-Term

- Treating physician and referring pulmonologist should decide together on optimal approach for follow-up
 - Recommended follow-up at 1, 3, 6 months and yearly
 - Evaluate patient for outcomes and complications
- Restart maintenance pulmonary rehabilitation at outpatient center and maintenance at home
- Treating physician should request that referring physician send patient back for re-evaluation if there is any of the following:
 - Loss of effect or no effect
 - Sudden loss of volume reduction on CT scan
 - Persistent cough
 - Persistent hemoptysis
 - Persistent or recurrent pneumonia
 - Frequent exacerbations not consistent with patients usual pre-EBV course

INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



Post-Procedural Re-Evaluation

Re-Evaluation Testing

- Recommend an additional low-dose HRCT scan to assess valve positioning if there is:
 - Breathing deterioration
 - No improvement
 - Cough
 - Hemoptysis
- Recommend the addition of bronchoscopic evaluation, valve adjustment or replacement if there is:
 - No or < 50% target lobe volume reduction on CT scan (at scheduled 30-45 day follow-up)
 - Sudden loss of benefit/loss of volume reduction on CT scan
 - Persistent cough
 - Persistent hemoptysis



INITIAL PATIENT
SCREENING

MEDICAL PATIENT
SELECTION

RADIOGRAPHIC
PATIENT SELECTION

PRE-ADMISSION
PROTOCOL

PROCEDURE DAY
PROTOCOL

POST-PROCEDURE
MANAGEMENT

PATIENT
FOLLOW-UP

POST-PROCEDURAL
RE-EVALUATION



References

- Criner GJ, Delage A, Voelker K. Late Breaking Abstract – Endobronchial Valves for Severe Emphysema – 12-month Results of the EMPROVE Trial. Eur Respir J. 2018;52(suppl 62). doi:10.1183/13993003.congress-2018.OA4928
- Herth FJF, Slebos D-J, Criner GJ, Shah PL. Endoscopic Lung Volume Reduction: An Expert Panel Recommendation - Update 2017. Respiration. 2017; 94(4):380-388. Doi: 10.1159/000479379
- Criner GJ. Disease and Anatomy Medical Management. Oral presentation at: Spiration Valve System - Treatment of Severe Emphysema Professional Education Course; June, 2019 (ongoing); Chicago, IL.
- Majid A. Management of Complications. Oral presentation at: Spiration Valve System - Treatment of Severe Emphysema Professional Education Course; June, 2019 (ongoing); Chicago, IL.
- Majid A. Patient Path. Oral presentation at: Spiration Valve System - Treatment of Severe Emphysema Professional Education Course; June, 2019 (ongoing); Chicago, IL.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Olympus is a registered trademark of Olympus Corporation, Olympus America Inc., and/or their affiliates. | Medical devices listed may not be available for sale in all countries.



OLYMPUS AMERICA INC.
3500 Corporate Parkway, PO Box 610, Center Valley, PA 18034

For more information, contact your Olympus sales representative, or call 800-848-9024.
www.medical.olympusamerica.com

©2019 Olympus America Inc. All rights reserved.
Printed in USA OAIRES0419EL30592