

## Microsoft® Update and PHP Update Monthly Notification

**Summary:** This notification contains information about Microsoft security updates released this month and approved for implementation on the Image Stream Medical EasySuite®, LiveStream®, nCare®/nStream®, and VaultStream® platforms. The only updates listed are those considered to be of a security or critical nature. Only the updates listed below have been validated with the EasySuite an EasySuite 4K, nCare GX/MX/RX, nStream G3, EasyLink Router, Thin Client, LiveStream StatusBoard Engine, MedPresence Stand-Alone, All-In-One Client, and VaultStream systems.

Also included is information on the latest PHP updates that are approved for the v8.x and up VaultStream systems.

**Products:** EasySuites, nStreams, nCares, and Thin Clients connected to a network and VaultStream, as well as the LiveStream StatusBoard Engine.

**Date:** December 10, 2019

### **Critical Security Bulletins (Obtained from Microsoft):**

**NOTE:** Image Stream Medical **does not approve** MS updates **KB975560, KB975562, KB2631813 or KB2780091** and strongly suggests they be excluded from automatic update scenarios. Failure to exclude these updates will result in equipment performance issues.

### **Approved Updates:**

<b>November and December 2019 Updates</b>		
<b>KB Number</b>	<b>Platform</b>	<b>Title</b>
<b>nStream v5.4 Updates</b>		
KB4525235	nStream	2019-11 Security Monthly Quality Rollup for Windows 7 for x86-based Systems
KB4523206	nStream	2019-11 Servicing Stack Update for Windows 7 for x86-based Systems
KB890830	nStream	Windows Malicious Software Removal Tool - December 2019
KB4530734	nStream	2019-12 Security Monthly Quality Rollup for Windows 7 for x86-based Systems
KB4531786	nStream	2019-12 Servicing Stack Update for Windows 7 for x86-based Systems
<b>VaultStream v5.3, 8.0.x, v9.0.x Updates – Windows® Server 2008</b>		
KB4523206	VaultStream 2008	2019-11 Servicing Stack Update for Windows Server 2008 R2 for x64-based Systems
KB4525235	VaultStream 2008	2019-11 Security Monthly Quality Rollup for Windows Server 2008 R2 for x64-based Systems
KB4520406	VaultStream 2008	2019-10 Preview of Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 7 and Server 2008 R2 for x64
KB890830	VaultStream 2008	Windows Malicious Software Removal Tool x64 - November 2019
KB890830	VaultStream 2008	Windows Malicious Software Removal Tool x64 - December 2019
KB4533095	VaultStream 2008	2019-12 Security and Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 7 and Server 2008 R2 for x64
KB4530734	VaultStream 2008	2019-12 Security Monthly Quality Rollup for Windows Server 2008 R2 for x64-based Systems
KB4531786	VaultStream 2008	2019-12 Servicing Stack Update for Windows Server 2008 R2 for x64-based Systems
<b>VaultStream v10.x Updates – Windows® Server 2012</b>		
KB4524445	VaultStream 2012	2019-11 Servicing Stack Update for Windows Server 2012 R2 for x64-based Systems
KB4525243	VaultStream 2012	2019-11 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems

KB890830	VaultStream 2012	Windows Malicious Software Removal Tool x64 - November 2019
KB4520408	VaultStream 2012	2019-10 Preview of Quality Rollup for .NET Framework 3.5, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 8.1 and Server 2012 R2 for x64
KB4520012	VaultStream 2012	2019-10 Preview of Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
B4533097	VaultStream 2012	2019-12 Security and Quality Rollup for .NET Framework 3.5, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 8.1 and Server 2012 R2 for x64
KB4530702	VaultStream 2012	2019-12 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
KB890830	VaultStream 2012	Windows Malicious Software Removal Tool x64 - December 2019
<b>VaultStream v10.x Updates – Windows® Server 2016</b>		
KB890830	VaultStream 2016	Windows Malicious Software Removal Tool x64 - November 2019
KB4525236	VaultStream 2016	2019-11 Cumulative Update for Windows Server 2016 for x64-based Systems
KB4520724	VaultStream 2016	2019-11 Servicing Stack Update for Windows Server 2016 for x64-based Systems
KB4519979	VaultStream 2016	2019-10 Cumulative Update for Windows Server 2016 for x64-based Systems
	VaultStream 2016	VMware, Inc. - System - 7/11/2019 12:00:00 AM - 9.8.16.0
KB890830	VaultStream 2016	Windows Malicious Software Removal Tool x64 - December 2019
KB4530689	VaultStream 2016	2019-12 Cumulative Update for Windows Server 2016 for x64-based Systems
<b>EasyLink Router v8.0 Updates</b>		
None		
<b>nCare RX v8.0 and v8.1, nCare MX/ nStream GX v8.0, v9.0, and Status Board Engine, 8.0 and v8.1 Updates</b>		
KB4519976	nCare RX/MX-SBE	2019-10 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x64-based Systems
KB4523206	nCare RX/MX-SBE	2019-11 Servicing Stack Update for Windows Embedded Standard 7 for x64-based Systems
KB4525235	nCare RX/MX-SBE	2019-11 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x64-based Systems
KB4525233	nCare RX/MX-SBE	2019-11 Security Only Quality Update for Windows Embedded Standard 7 for x64-based Systems
KB4524741	nCare RX/MX-SBE	2019-11 Preview of Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 7 and Server 2008 R2 for x64
		windows6.1-kb4019990-x64_35cc310e81ef23439ba0ec1f11d7b71dd34adfe5.msu
		windows6.1-kb4507004-x64_a7d988d65422c8cfbbbed141c95aeea37a8743167.msu
KB4524421	nCare RX/MX-SBE	2019-11 Preview of Quality Rollup for .NET Framework 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64
KB4531182	nCare RX/MX-SBE	2019-11 Preview of Quality Rollup for .NET Framework 4.8 for Windows 7 and Server 2008 R2 for x64
KB4533095	nCare RX/MX-SBE	2019-12 Security and Quality Rollup for .NET Framework 3.5.1, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows 7 and Server 2008 R2 for x64
		windows6.1-kb4019990-x64_35cc310e81ef23439ba0ec1f11d7b71dd34adfe5.msu



**Image Stream Medical**

One Monarch Drive Littleton, MA 01460 U.S.A.  
[www.olympusamerica.com/olympus-systems-integration](http://www.olympusamerica.com/olympus-systems-integration)

Image Stream Medical, Inc., is a wholly owned subsidiary of Olympus Corporation of the Americas

		windows6.1-kb4507004-x64_a7d988d65422c8cfbbbed141c95aeea37a8743167.msu
KB4533012	nCare RX/MX-SBE	2019-12 Security and Quality Rollup for .NET Framework 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2 for Windows 7 and Server 2008 R2 for x64
KB4533005	nCare RX/MX-SBE	2019-12 Security and Quality Rollup for .NET Framework 4.8 for Windows 7 and Server 2008 R2 for x64
KB4531786	nCare RX/MX-SBE	2019-12 Servicing Stack Update for Windows Embedded Standard 7 for x64-based Systems
KB4530734	nCare RX/MX-SBE	2019-12 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x64-based Systems
KB4530692	nCare RX/MX-SBE	2019-12 Security Only Quality Update for Windows Embedded Standard 7 for x64-based Systems
<b>HP Thin Client v8.0 and EasyLink Router v8.1 Updates</b>		
KB4519976	Thin Client v8.0-ELR v8.1	2019-10 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x86-based Systems
KB4523206	Thin Client v8.0-ELR v8.1	2019-11 Servicing Stack Update for Windows Embedded Standard 7 for x86-based Systems
KB4525235	Thin Client v8.0-ELR v8.1	2019-11 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x86-based Systems
KB4525233	Thin Client v8.0-ELR v8.1	2019-11 Security Only Quality Update for Windows Embedded Standard 7 for x86-based Systems
KB4531786	Thin Client v8.0-ELR v8.1	2019-12 Servicing Stack Update for Windows Embedded Standard 7 for x86-based Systems
KB4530734	Thin Client v8.0-ELR v8.1	2019-12 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x86-based Systems
KB4530692	Thin Client v8.0-ELR v8.1	2019-12 Security Only Quality Update for Windows Embedded Standard 7 for x86-based Systems
<b>nCare MX v10.2.1</b>		
KB4520724	nCare MX	2019-11 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB4525236	nCare MX	2019-11 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4530689	nCare MX	2019-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
<b>MedPresence Stand-Alone v9.0.x and All-In-One Client (works with v9.0.x and 10.x)</b>		
KB4516115	MedPresence SA, and All-In-One Client	2019-09 Security Update for Adobe Flash Player for Windows 10 Version 1607
<b>EasySuite 4K v10.1.2   EasySuite 4K Controller (a.k.a. BlizCon)</b>		
KB4525236	EasySuite 4K Controller	2019-11 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4520724	EasySuite 4K Controller	2019-11 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB4530689	EasySuite 4K Controller	2019-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
<b>EasySuite 4K v10.1.2   EasySuite 4K Video Manager (a.k.a. nCare RX2)</b>		
KB4519998	EasySuite 4K Video Manager	2019-10 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4521858	EasySuite 4K Video Manager	2019-10 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB4525236	EasySuite 4K Video Manager	2019-11 Cumulative Update for Windows 10 Version 1607 for x64-based Systems



**Image Stream Medical**

One Monarch Drive Littleton, MA 01460 U.S.A.  
[www.olympusamerica.com/olympus-systems-integration](http://www.olympusamerica.com/olympus-systems-integration)

Image Stream Medical, Inc., is a wholly owned subsidiary of Olympus Corporation of the Americas

KB4520724	EasySuite 4K Video Manager	2019-11 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB4530689	EasySuite 4K Video Manager	2019-12 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
<b>EasySuite 4K v10.0.x – LiveStream StatusBoard 2</b>		
KB890830	LiveStream StatusBoard 2	Windows Malicious Software Removal Tool x64 - November 2019
KB4525236	LiveStream StatusBoard 2	2019-11 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4520724	LiveStream StatusBoard 2	2019-11 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB4519979	LiveStream StatusBoard 2	2019-10 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB4519562	LiveStream StatusBoard 2	2019-10 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB4521858	LiveStream StatusBoard 2	2019-10 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB4524095	LiveStream StatusBoard 2	2019-10 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
KB890830	LiveStream StatusBoard 2	Windows Malicious Software Removal Tool x64 - October 2019
<b>PHP Updates</b>		
5.6.40 and 7.1.33		

For further questions concerning this notice, please call Technical Support at (866) 376-6800



**Image Stream Medical**

One Monarch Drive Littleton, MA 01460 U.S.A.  
[www.olympusamerica.com/olympus-systems-integration](http://www.olympusamerica.com/olympus-systems-integration)

Image Stream Medical, Inc., is a wholly owned subsidiary of Olympus Corporation of the Americas