

Microsoft® Update and PHP Update Monthly Notification

Summary: This notification contains information about Microsoft security updates released this month and approved for implementation on the Gyrus ACMI, Inc., ES4K, LiveStream™, nCare™, and VaultStream™ platforms. The only updates listed are those considered to be of a security or critical nature. Only the updates listed below have been validated with the EasySuite HD and ES4K, nCare MX/RX, EasyLink Router, LiveStream StatusBoard Engine, and VaultStream systems.

Also included is information on the latest PHP updates that are approved for the v10.x and up VaultStream systems.

Products: EasySuite HD, ES4K and nCares, connected to a network and VaultStream, as well as the LiveStream StatusBoard Engine.

Date: September 14, 2021

Critical Security Bulletins (Obtained from Microsoft):

NOTE: Gyrus ACMI, Inc. **does not approve** MS updates **KB975560, KB975562, KB2631813 or KB2780091** and strongly suggests they be excluded from automatic update scenarios. Failure to exclude these updates will result in equipment performance issues.

Approved Updates:

September 2021 Updates		
KB Number	Platform	Title
VaultStream v10.x Updates – Windows® Server 2012		
KB5005613	VaultStream 2012	2021-09 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems
KB890830	VaultStream 2012	Windows Malicious Software Removal Tool x64 - v5.93
VaultStream v10.x Updates – Windows® Server 2016		
KB890830	VaultStream 2016	Windows Malicious Software Removal Tool x64 - v5.93
KB5005573	VaultStream 2016	2021-09 Cumulative Update for Windows Server 2016 for x64-based Systems
KB5005698	VaultStream 2016	2021-09 Servicing Stack Update for Windows Server 2016 for x64-based Systems
VaultStream v10.4 Updates – Windows® Server 2019		
KB5005568	VaultStream 2019	2021-09 Cumulative Update for Windows Server 2019 (1809) for x64-based Systems
KB5005102	VaultStream 2019	2021-08 Cumulative update Preview for Windows Server 2019 (1809) for x64-based Systems
KB890830	VaultStream 2019	Windows Malicious Software Removal Tool x64 - v5.93
KB2267602	VaultStream 2019	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.349.801.0)
KB2267602	VaultStream 2019	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.349.778.0)
KB2267602	VaultStream 2019	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.349.720.0)
nCare MX v10.x		
KB5005698	nCare MX	2021-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB5005573	nCare MX	2021-09 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB5005043	nCare MX	2021-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems

EasySuite 4K v10.x E4K Controller (a.k.a. BlizCon)		
KB5005698	ES4K Controller	2021-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB5005573	ES4K Controller	2021-09 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
EasySuite 4K v10.x ES4K Video Manager (a.k.a. nCare RX2)		
KB5005698	ES4K Video Manager	2021-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB5005573	ES4K Video Manager	2021-09 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
EasySuite 4K v10.x – LiveStream StatusBoard 2		
KB5005112	LiveStream StatusBoard 2	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems
KB2267602	LiveStream StatusBoard 2	Security Intelligence Update for Microsoft Defender Antivirus (Version 1.349.811.0)
MedPresence Tower v10.3		
KB5005698	MedPresence Tower	2021-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems
KB5005573	MedPresence Tower	2021-09 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB5005043	MedPresence Tower	2021-08 Cumulative Update for Windows 10 Version 1607 for x64-based Systems
KB5004752	MedPresence Tower	2021-08 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64
PHP Updates		
7.1.33 and 7.4.16 – shipping versions.		
7.3.30 is available upon request for v10.2.2 VaultStreams only.		
7.4.23 is available upon request for v10.4 VaultStreams only.		

For further questions concerning this notice, please call Technical Support at 1-800-848-9024